



Hard on Pathogens, Safe on Surfaces.

A powerful, one-step cleaner/disinfectant effective against a broad spectrum of pathogens, including HIV, H1N1, MRSA and many others. Useful in disinfecting vinyl exam tables, counter tops and hard surfaces, mammography compressor plates, ultrasound transducers and other hard, nonporous, non-critical surfaces.

PROTEX® DISINFECTANTS CAN BE USED AGAINST 2019 NOVEL CORONAVIRUS

Protex® Disinfectant Spray demonstrated effectiveness against viruses similar to the 2019 Novel Coronavirus on hard, non-porous surfaces. Therefore Protex can be used against the 2019 Novel Coronavirus when used in accordance with the Directions for Use.

FEATURES



Effective

Protex® disinfectant spray is effective against more than 30 different pathogens including HIV, H1N1, MRSA and more.

Safe on Surfaces

Alcohol, bleach and solvent free disinfectant is gentle on sensitive equipment.

Convenient

32oz trigger spray bottle for use in both hospital and institutional areas.

PROTEX® DISINFECTANT SPRAY

Product #	Cross Ref	Description	# / Case
32021041	42-32	32 oz trigger spray bottle	12



proteXTM disinfectant
spray

Active Ingredients:

Octyl decyl dimethyl ammonium chloride..... 0.025%

Dioctyl dimethyl ammonium chloride 0.010%

Didecyl dimethyl ammonium chloride 0.015%

Alkyl (C₁₄, 50%; C₁₂, 40%; C₁₆, 10%)

dimethyl benzyl ammonium chloride..... 0.034%

Other Ingredients 99.916%

Total..... 100.000%

LA-032-0000F

KEEP OUT OF REACH OF CHILDREN

CAUTION

See inside panel for Precautionary Statements and
First Aid


EPA reg. No. 6836-152-82613

EPA Establishment No. 82613-NJ-001

Parker Laboratories Inc.
286 Eldridge Road
Fairfield, NJ 07004

www.parkerlabs.com

NET CONTENTS 32 FL OZ



proteXTM disinfectant
spray

REF 42-32

One-Step
Disinfectant

For Hospitals and Institutional Use

Effective in the presence of
5% serum contamination

- Kills Pandemic 2009 H1N1 Influenza A Virus
- Staphylocidal
- Pseudomonacidal
- Bactericidal
- Salmonellacidal
- Fungicidal
- Virucidal*



PEEL
HERE

USE PROTEX™ DISINFECTANT SPRAY

- As a disinfectant on hard, nonporous surfaces
- For disinfecting of ultrasound transducers, probes, mammography compressor plates and other hard nonporous surfaces
- To clean and disinfect tables, chairs and countertops
- Use on nonporous athletic mats, wrestling mats, gymnastic mats, exercise equipment and training tables

USE IN:

- Hospitals • Medical Related Facilities • Nursing Homes
- Athletic Facilities • Exercise Facilities

USE ON:

- Washable, hard, nonporous surfaces of: walls, counters, ultrasonic baths
- Other hard nonporous surfaces made of: laminated surfaces, metal, plastic, stainless steel, vinyl and plastic upholstery, finished woodwork

Protex Disinfectant Spray is a one-step hospital use disinfectant, Bactericidal according to the current AOAC Use-Dilution Test Method, Fungicidal according to the AOAC Fungicidal Test and Virucidal* according to the Virucidal qualification modified in the presence of 5% organic serum against:

Bacteria: *Acinetobacter baumannii* • *Corynebacterium ammoniagenes* • *Burkholderia cepacia* • *Campylobacter jejuni* • *Enterobacter aerogenes* • *Enterobacter cloacae* NDM-1 • *Enterococcus faecalis*—Vancomycin resistant—VRE • ESBL *Escherichia coli* • *Escherichia coli*, New Delhi Metallo-Beta Lactamase (NDM-1) • *Klebsiella pneumoniae* • *Klebsiella*

pneumoniae—NDM-1 positive • *Legionella pneumophila* • *Pseudomonas aeruginosa* • *Salmonella enterica* • *Salmonella schottmuelleri* • *Salmonella typhi* • *Serratia marcescens* • *Shigella dysenteriae* • *Staphylococcus aureus*—Community Associated Methicillin-Resistant [CA-MRSA] [NRS384][USA300] • *Staphylococcus aureus*—Methicillin Resistant—MRSA • *Staphylococcus aureus* Multi-Drug Resistant • *Staphylococcus aureus*—Vancomycin Intermediate Resistant—VISA • *Streptococcus pyogenes* • *Vibrio cholerae*

Viruses: *Hepatitis B Virus (HBV) • *Hepatitis C Virus (HCV) • *Herpes Simplex Virus Type 1 • *Herpes Simplex Virus Type 2 • *HIV-1 (AIDS Virus) • *Human Coronavirus • *Influenza Type A / Brazil • *2013 Influenza A Virus (H7N9) • *Norwalk Virus • *Respiratory Syncytial Virus • *Rotavirus • *SARS Associated Coronavirus • *Vaccinia

Fungi: *Trichophyton mentagrophytes* • *Aspergillus niger* • *Aspergillus fumigatus*

DIRECTIONS FOR USE

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

DISINFECTION / *VIRUCIDAL / FUNGICIDAL DIRECTIONS:

Spray 6-8 inches from surface, until surfaces are thoroughly wet. Do not breathe spray. Treated surfaces must remain wet for 10 minutes. For Influenza Virus Type A and Human Coronavirus, treated surfaces must remain wet for 1 minute. Wipe dry with a clean cloth, sponge, or mop or allow to air dry. For heavily soiled areas, a pre-cleaning is required.

Rinse all surfaces that come in contact with food such as countertops, appliances, tables and stovetops with potable water before reuse. Do not use on utensils, glassware and dishes.

For fungicidal activity: Protex Disinfectant Spray is an effective one-step fungicide against *Trichophyton mentagrophytes* when used on hard nonporous surfaces in areas such as locker rooms, dressing rooms, shower and bath areas and exercise facilities. Allow surface to remain wet for 10 minutes.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body.

ULTRASONIC BATH DISINFECTANT / *VIRUCIDAL

DIRECTIONS: Use Protex Disinfectant Spray to disinfect hard nonporous non-critical objects compatible with Ultrasonic cleaning units. Pour Protex Disinfectant Spray directly into bath chamber. Pre-clean heavily soiled objects. Place objects into unit and operate for a minimum of 10 minutes, according to the manufacturer's directions. Remove objects and wash with sterile water (sterile water for injection) or allow to air dry. Use fresh product for each use.

KILLS HIV-1 AND HBV AND HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES / OBJECTS PREVIOUSLY SOILED WITH BLOOD / BODY FLUIDS

in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces / objects with blood or body fluids and in which the surfaces / objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS) or Hepatitis B Virus (HBV) or Hepatitis C Virus (HCV).

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1 OR HBV OR HCV ON SURFACES / OBJECTS SOILED WITH BLOOD / BODY FLUIDS.

Personal protection: Clean-up should always be done wearing protective gloves, gowns, masks, and eye protection.

Cleaning Procedure: Blood and other body fluids containing HIV or HBV or HCV must be thoroughly cleaned from surfaces and objects before application of Protex Disinfectant Spray.

Contact time: Leave surface wet for 1 minute (60 seconds) for HIV-1 and 10 minutes for HBV and HCV. Use a 10 minute contact time for disinfection against all other viruses, bacteria and fungi claimed.

Disposal of Infectious material: Blood, body fluids, cleaning materials and clothing should be autoclaved and disposed of according to local regulations for infectious waste disposal.

PROTEX™ DISINFECTANT SPRAY

Is a one-step Hospital Use disinfectant cleaner that is effective against a broad spectrum of bacteria, is virucidal* including HIV-1, the AIDS Virus and Hepatitis B Virus and Hepatitis C Virus, and fungicidal when used as directed.

- Meets AOAC efficacy standards for hospital Disinfectants

USE PROTEX DISINFECTANT SPRAY

- As a disinfectant on hard, nonporous surfaces
- For disinfecting of ultrasound transducers, probes, mammography compressor plates and other hard nonporous surfaces

PRECAUTIONARY STATEMENTS

HAZARDS TO HUMANS AND DOMESTIC ANIMALS

Caution. Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, or using tobacco.

Parker Laboratories Inc.
286 Eldridge Road
Fairfield, NJ 07004

www.parkerlabs.com

FIRST AID

IF IN EYES: Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present after first 5 minutes, then continue rinsing the eye.

Call a poison control center or doctor for treatment advice. Have the container or label with you when calling a poison control center or doctor, or going for treatment.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal.

Pesticide Storage:

Open dumping is prohibited. Store in original container in areas inaccessible to children.

Container Disposal:

Nonrefillable container. Do not reuse or refill this container. Triple rinse (or equivalent). Offer for recycling or reconditioning or puncture and dispose of in a sanitary landfill, or incineration, or if allowed by state and local authorities, by burning. If burned, stay out of smoke.

NET CONTENTS 32 FL OZ

disinfectant
spray
protex™

REF 42-32

One-Step Disinfectant

For Hospitals and Institutional Use

Effective in the presence of
5% serum contamination

- Kills Pandemic 2009 H1N1 Influenza A Virus
- Staphylocidal
- Pseudomonacidal
- Bactericidal
- Salmonellacidal
- Fungicidal
- Virucidal*



SECTION 1 — CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

Product Name: PROTEX® DISINFECTANT SPRAY
(As used on label and List)

Product Use: BIOCIDES

Manufacturer's Name: PARKER LABORATORIES INC.
Address: 286 ELDRIDGE ROAD
(number, street, City, State, Zip code) FAIRFIELD, NJ 07004 USA

Emergency Telephone Number: (800)-631-8888
Telephone Number for Information: (973)-276-9500 (8:30am – 5:15pm ET)

SECTION 2 — HAZARDS IDENTIFICATION

Hazardous Components (Specific Chemical Identity; Common Name(s)): NONE

Classification of the Substance/Mixture: NOT CLASSIFIED

Signal Word: N/A

Hazard Statements: NO KNOWN SIGNIFICANT EFFECTS OR CRITICAL HAZARDS

Precautionary Statement: HANDLE IN ACCORDANCE WITH GOOD INDUSTRIAL HYGIENE PRACTICES

OSHA Hazards: NO OSHA HAZARDS

SECTION 3 — COMPOSITION/INFORMATION ON INGREDIENTS

Hazardous Ingredients (specific)	% Composition	CAS Number
N/A	N/A	N/A

SECTION 4 — FIRST AID MEASURES

Signs and Symptoms of Exposure: N/A

Medical Symptoms Generally Aggravated by Exposure: NONE KNOWN

Emergency and First Aid Procedures:

Inhalation: Move to fresh air immediately. If experiencing difficulty breathing, seek medical attention.

Skin contact: In the case of skin irritation or allergic reactions, wash with soap and plenty of water. See a physician.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Ingestion: If accidentally swallowed obtain immediate medical attention. Clean mouth with water and drink afterwards plenty of water.

SAFETY DATA SHEET

Product Name: PROTEX® DISINFECTANT SPRAY

Do not induce vomiting without medical advice.
Never give anything by mouth to an unconscious person.

SECTION 5 — FIRE FIGHTING MEASURES

Flash Point(Method Used): >199.9 °F **Flammable Limits:** N/A **LEL:** N/A **UEL:** N/A

Extinguishing Media: DRY POWDER, WATER SPRAY, FOAM. USE WATER SPRAY TO COOL UNOPENED CONTAINERS

Special Fire Fighting Procedures: NONE

Unusual Fire and Explosion Hazards: HEATING OR FIRE CAN RELEASE TOXIC GAS

Protective Equipment and Precautions for Firefighters: WEAR A SELF-CONTAINED BREATHING APPARATUS AND FULL PROTECTIVE GEAR.

SECTION 6— ACCIDENTAL RELEASE MEASURES

Leak and Spill Procedures: SWEEP UP OR VACUUM UP SPILLAGE AND COLLECT IN SUITABLE CONTAINER FOR DISPOSAL.

Personal Precautions: USE RESPIRATOR WHEN PERFORMING OPERATIONS INVOLVING POTENTIAL EXPOSURE TO VAPOR OF THE PRODUCT

Environmental Precautions: PREVENT PRODUCT FROM ENTERING DRAINS

SECTION 7— HANDLING AND STORAGE

Handling procedures and Equipment: PROVIDE SUFFICIENT AIR EXCHANGE AND/OR EXHAUST IN WORK ROOMS.

Storage requirements: KEEP CONTAINER TIGHTLY CLOSED AWAY FROM FOOD, DRINKS AND ANIMAL FEEDING STUFF. STORAGE TEMPERATURE <60°C

SECTION 8— EXPOSURE CONTROL/PERSONAL PROTECTION

Exposure Limits: ☐ ACGIH TLV ☐ OSHA PEL ☐ Other(Specify)

Specific Engineering Controls (such as ventilation, enclosed process) ENSURE ADEQUATE VENTILATION

Personal Protective Equipment:

☒ Gloves ☒ Eye ☐ Clothing ☐ Footwear ☐ Respirator ☐ Other (specify)

If Checked, please specify type: HAND PROTECTION NOT REQUIRED WHEN USED AS DIRECTED, AVOID EYE CONTACT. EMERGENCY RESPONDERS SHOULD WEAR TIGHT FITTING SAFETY GOGGLES AND IMPERMEABLE GLOVES



Physical State(Appearance):	Clear Liquid
Color:	Blue Tint
Odor:	Mild Lemon
Specific Gravity(H₂O=1):	N/D
pH:	10.0 - 12.0
Solubility in Water:	Soluble
Boiling point:	N/D
Melting Point:	N/D
Flammability(solid, gas):	N/D
Vapor Pressure(mm Hg):	N/D
Vapor Density(AIR=1):	N/D
Viscosity:	N/A
Density:	N/D
Evaporation Rate(Butyl Acetate= 1):	N/D

Stability: Unstable ☐ Conditions to Avoid: NONE
Stable ☒

Incompatibility (Materials to Avoid): NONE

Conditions to Avoid: NONE

Reactivity: NO SPECIFIC TEST DATA RELATED TO REACTIVITY AVAILABLE FOR THIS PRODUCT.

Health Hazards (Acute and Chronic): NONE KNOWN

Route(s) of Entry: ☐ Skin Contact ☐ Eye Contact ☐ Inhalation ☒ Ingestion

Carcinogenicity: NTP? NO IARC Monographs? NO OSHA Regulated: NO

Signs and Symptoms of Exposure: N/A

Medical Conditions Generally Aggravated by Exposure: NONE KNOWN

SECTION 12 — ECOLOGICAL INFORMATION

Ecotoxicity:	N/A
Persistence and Degradability	THERE IS NO DATA AVAILABLE FOR THIS PRODUCT
Bioaccumulative Potential	THERE IS NO DATA AVAILABLE FOR THIS PRODUCT
Mobility in the Soil:	THERE IS NO DATA AVAILABLE FOR THIS PRODUCT
Additional Ecological Information:	THERE IS NO DATA AVAILABLE FOR THIS PRODUCT

SECTION 13 — DISPOSAL CONSIDERATIONS**Waste Disposal:**

FOLLOW ALL FEDERAL, STATE AND LOCAL REGULATIONS FOR NON-HAZARDOUS WASTE DISPOSAL.

SECTION 14 — TRANSPORT INFORMATION

DOT Classification:	NOT REGULATED
IMDG Classification:	NOT REGULATED
TDG Classification:	NOT REGULATED
IATA Classification:	NOT REGULATED
RID:	NOT REGULATED

SECTION 15 — REGULATORY INFORMATION

OSHA Hazards:	NO OSHA HAZARDS
SARA 311/312 Hazards:	NO SARA HAZARDS
California Prop. 65:	This product does not contain any chemicals known to the State of California to cause cancer, birth defects or other reproductive harm.

EPA Registration Number: 6836-152-82613

SECTION 16 — OTHER INFORMATION

Date of Issue:	May 15 2020
Prepared By:	Quality Assurance
Version:	2

Disclaimer: This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. No warranty, guarantee or representation is made as to the correctness and sufficiency of the information. It's the user's responsibility to determine safe conditions for the use of this product and to assume liability for loss, injury, or expense resulting from the misuse of this product.